

Rabbit Gamma Globulin

Product Description

Rabbit gamma globulin is purified from normal rabbit serum (either sex, healthy animals) using different fractionation procedures and lyophilisation.

Applications

Useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.

Physical State

Off white fine, flaky powder

Testing

Serum Protein Electrophoresis: $\geq 95\%$

Protein by Weight: O.D. 280 nm E1%1.39: $\geq 85\%$ - Residual salts, lipids, etc., comprise the rest of the weight

Moisture Content: Report result

Bioburden: ≤ 1000 CFU/g

Packaging

10 g: Polypropylene bottle; 500 g: HDPE container with lid

Storage Temperature

+ 2-8°C

Shipping Conditions

Ambient temperature

Applications

Useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.

Product Code

Code	Description	Size
GG-005-10g	Rabbit Gamma Globulin - US origin	10g
GG-005-100g	Rabbit Gamma Globulin - US origin	100g
GG-005-500g	Rabbit Gamma Globulin - US origin	500g

For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Group is [ISIA Traceability Certified](#)



To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk