

Sheep Gamma Globulin

Product Description

Sheep gamma globulin is purified from normal sheep serum or plasma using a combination of proprietary methods.

Physical State

Off white fine, flaky powder

Testing

Serum Protein Electrophoresis: $\geq 98\%$

Protein by Weight: O.D. 280 nm $E_{1\%}^{1.39}$: $\geq 85\%$ - Residual sales, lipids, etc comprise the rest of the weight

Moisture Content: FIO

Packaging

Polypropylene bottle

Storage Temperature

+ 2-8°C

Shipping Conditions

Ambient temperature

Applications

Useful as a stabilizer/blocker for EIA and RIA systems, IHC assays, immunoprecipitation, Western Blot and ELISA experiments, and protein standards.

Product Code

Code	Description	Size
GG-008-1g	Sheep Gamma Globulin – US origin	1g
GG-008-10g	Sheep Gamma Globulin - US origin	10g
GG-008-50g	Sheep Gamma Globulin - US origin	50g

For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Group is [ISIA Traceability Certified](#)



To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk