

# Foetal Bovine Serum

LSP offer South American FBS from Brazil (also referred to as EU Grade FBS). This material is collected under specific protocols and imported into the EU under Commission Regulation (EU) No 142/2011. The material is processed in the UK at LSP facilities, using fully disposable systems.

## FEATURES

Our quality assurance process means each batch of FBS is fully tested prior to release. Full QA Certification supplied.

- **Filtered:** triple filtered to 0.1 µm
- **Sterile tested:** for the absence of bacteria, fungi, yeast and mycoplasma
- **Viral tested:** for bovine diseases BVD-V, PI3, IBR and BVB-AB
- **Biological performance:** is assessed for cell growth, plating efficiency and cloning efficiency
- **Endotoxin:** levels are assessed
- **Traceability certified:** LSP is a member of the International Serum Industry Association (ISIA) and is ISIA Traceability certified.



## BENEFITS

- Competitively priced
- A range of batches are available to test from
- Batch reservations are available FOC - sample size is 50mL
- No minimum order quantities



### Ordering Information

S-001A-BR	Foetal Bovine Serum, South American origin, Aseptically filled, 500mL
S-001B-BR	Foetal Bovine Serum, South American origin, Aseptically filled, 100mL
S-001A-HI-BR	Foetal Bovine Serum, South American origin, Aseptically filled, Heat Inactivated, 500mL
S-001B-HI-BR	Foetal Bovine Serum, South American origin, Aseptically filled Heat Inactivated, 100mL

Please note there is a shipping /dry ice charge per delivery.

#### To request a quote or to batch test:

T: 01234 889180 or email: [sales@lifescienceproduction.co.uk](mailto:sales@lifescienceproduction.co.uk)

Web: [www.lifescienceproduction.co.uk/products/](http://www.lifescienceproduction.co.uk/products/)

Follow us: @LifeScienceGrp

Life Science Production is a division of Life Science Group Ltd. Life Science Group Ltd is an ISO 9001:2015 and ISIA Traceability Certified company.

**LSP**  
LIFE SCIENCE PRODUCTION