

BSA Polymerised Type IV – lyophilised

Product Description

Bovine Serum Albumin (BSA) is a non-glycosylated protein of 66 kDa. Produced by the liver, it is the most abundant protein in plasma.

LSP Bovine Serum Albumin is obtained by a proprietary Heat-Shock method, designed to prevent the excessive use of denaturing organic solvents. It undergoes a thermal inactivation during 3 hours at 65 °C in order to inactivate potential viruses. This 30% solution BSA has been 0.2 µm filtered and is supplied as a ready-to-use solution in deionised water. This product contains 0.09% sodium azide.

Bovine plasma used for the manufacture of this product is collected from French healthy cattle declared free from BSE disease in its clinical forms. It is sourced in EU-approved French slaughterhouses.

Applications

Albumin has a key role for maintaining oncotic pressure, and is the main carrier for low molecular weight compounds such as fatty acids and hormones. Due to its versatile properties, BSA plays a major role in numerous biological applications such as blood grouping, cell culture or immunoassays (ELISA, RIA etc).

BSA Polymerised Type I solution has been specifically developed for use as a coating agent to avoid non specific interactions

Physical Characteristics

Appearance: Pale yellow flakes.

Analysis	Method	Specification
Dry extract	Loss on drying	> 300 g/L
Total proteins	Kjeldahl	> 290 g/L
pH	Potentiometry	6.5 - 7.4
Albumin purity	Electrophoresis	> 97 %
Mesophilic germs	Culture	Absence / 1 mL
Proteases	Protease activity assay	Undetectable
Monomers rate	Size exclusion chromatography	< 75 %

Product Information

Product Code	Product Description	Pack Size
BSA-7104470-1Kg	BSA Polymerized Type IV lyophilised - EU origin	1 Kg

Additional Information

CAS number: 9048-46-8

Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Production is [ISIA Traceability Certified](#)

Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk

Address: PO Box 1519, Bedford, United Kingdom

