

Antibiotic/Antimycotic 100x

Product Code: LSM0010B

Product Description

Inorganic salts

Sodium chloride 8500 mg/mL

Other components

Amphotericin B 25 µg/mL

Penicillin G 10,000 U/mL

Streptomycin sulphate 10 mg/mL

Presentation

100 mL

Storage

Store at -20°C

For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Production is [ISIA Traceability Certified](#)

Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk

Address: PO Box 1519, Bedford, United Kingdom