

Sodium Pyruvate Solution

Sodium Pyruvate solution is used in cell culture in the formulation of broths and media as a supplement to improve cell survival. It is used by cells as an easily accessible carbohydrate source. Additionally, it is involved with amino acid metabolism and initiates the Krebs's cycle.

Complete Certificate of Analysis available for each production lot, along with full formulation information.

Storage: 2°C to 8°C

Working dilution: 1:100 for most cell culture applications

Ordering information

Product Code	Product Description	Size
25-000-CI	Corning® 100 mL Sodium Pyruvate, Liquid 100 mM Solution [+] 8.5 g/L NaCl	6 x 100mL

For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Production is [ISIA Traceability Certified](#)

Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk

Address: PO Box 1519, Bedford, United Kingdom