

Gentamycin sulphate 10mg/mL

Product Code: LSM0011B/E

Product Description

Antibiotics, combined with good sterile technique, help to prevent microbiological contamination against Gram negative, Gram-positive bacteria and mycoplasma.

Gentamicin Sulphate is a broad spectrum bactericidal with activity against gram negative and gram positive organisms including Pseudomonas. It is stable over a wide range of pH and temperature and can be autoclaved. It's normal appearance may range from nearly colourless to amber or yellow.

Other components

Gentamycin sulphate	10000 mg/mL
Water	10,000 U/mL

Presentation:

LSM0011B	100 mL
LSM0011E	10 mL

Storage

Store at -20°C

For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.

Support

Life Science Production is a division of Life Science Group Ltd.
Life Science Production is [ISIA Traceability Certified](#)
Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: sales@lifesciencegroup.co.uk

Website: www.lifescienceproduction.co.uk

Address: PO Box 1519, Bedford, United Kingdom