

## Gentamycin sulphate 50mg/mL

Product Code: LSM0012E

### Product Description

Antibiotics, combined with good sterile technique, help to prevent microbiological contamination against Gram negative, Gram-positive bacteria and mycoplasma.

Gentamicin Sulphate is a broad spectrum bactericidal with activity against gram negative and gram positive organisms including Pseudomonas. It is stable over a wide range of pH and temperature and can be autoclaved. It's normal appearance may range from nearly colourless to amber or yellow.

#### Other components

Gentamycin sulphate	50000 mg/mL
Water	10,000 U/mL

### Presentation

LSM0012E 10 mL

### Storage

Store at -20°C

**For Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in clinical or *in vitro* procedures.**

### Support

Life Science Production is a division of Life Science Group Ltd.  
Life Science Production is [ISIA Traceability Certified](#)  
Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: [sales@lifesciencegroup.co.uk](mailto:sales@lifesciencegroup.co.uk)

Website: [www.lifescienceproduction.co.uk](http://www.lifescienceproduction.co.uk)

Address: PO Box 1519, Bedford, United Kingdom