

Corning® Transport Medium

An isotonic solution at neutral pH containing protective proteins, antibiotic to prevent microbial growth, and anti-mycotic to prevent contamination by yeast and fungi.

Formulated using the final ingredient concentrations outlined in CDC's SOP#: DSR-052-01: Preparation of Viral Transport Media.

Product Information

- Size: 100 mL each / 6 per case (current). Coming soon 500 mL & 1 L
- Shipping and storage: 2°C to 8°C
- Shelf life: 12 months
- Sterility: 10⁻³ SAL, labelled "Sterile A"
- Made in an FDA-registered (Class I Medical Devices) and ISO 13485:20013 certified facility
- cGMP grade

| Product Code | Product Description | Size |
|--------------|---------------------------|---------|
| 25-500-CI | Corning® Transport Medium | 100 mL |
| 25-500-CV | Corning® Transport Medium | 500 mL |
| 25-500-CM | Corning® Transport Medium | 1 litre |

Support

Life Science Production is a division of Life Science Group Ltd.
Life Science Production is [ISIA Traceability Certified](#)
Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:
Telephone: +44 (0) 1234 889180
Email: sales@lifesciencegroup.co.uk
Website: www.lifescienceproduction.co.uk
Address: PO Box 1519, Bedford, United Kingdom