

## LSG-1200B - L-Methionine 200mM

### Product Description

Methionine is an essential amino acid in humans. As the precursor of other amino acids such as cysteine and taurine, versatile compounds such as SAM-e, and the important antioxidant glutathione, methionine plays a critical role in the metabolism and health of many species, including humans.

### Product Specification

Components	Molecular Weight	Concentration
L- Methionine	149.2	200mM

### Product information

Product Code	Product Description	Size
LSG-1200B	L-Methionine 200mM	100mL

### Storage & Handling

Store frozen at -20°C or below

### Shipping

Ships on dry ice.

### Support

Life Science Production is a division of Life Science Group Ltd.

Life Science Production is [ISIA Traceability Certified](#)

Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:

Telephone: +44 (0) 1234 889180

Email: [sales@lifesciencegroup.co.uk](mailto:sales@lifesciencegroup.co.uk)

Website: [www.lifescienceproduction.co.uk](http://www.lifescienceproduction.co.uk)

Address: PO Box 1519, Bedford, United Kingdom



## QUALITY - TRACEABILITY - INTEGRITY

### Life Science Production

The Science Centre, 1 Blaydon Road, Sandy, Bedfordshire, SG19 1RZ

Tel: +44 (0) 1234 889180 Email: [sales@lifesciencegroup.co.uk](mailto:sales@lifesciencegroup.co.uk) Web: [www.lifesciencegroup.co.uk](http://www.lifesciencegroup.co.uk)

Life Science Production, Antibody Production Services and Diagnostic Kit Services are divisions of Life Science Group Ltd

Registered in England and Wales Company Number: 06606016

Page 1 of 1