

## Veronal Buffer

Veronal Buffer is used as a diluent and/or stabilizer in serological testing.

### Product Specification

Product description	Concentration
CaCl <sub>2</sub> ·2H <sub>2</sub> O	185 mg/L
MgCl <sub>2</sub> ·6H <sub>2</sub> O	840 mg/L
NaCl	42,500 mg/L
Sodium Barbitol	1,870 mg/L
Barbitol	2,870 mg/L

### Product information

Product Code	Product Description	Size
LSG-1174B	Veronal Buffer 5x	100mL

### Storage

Product is stored at between 2°C and 8°C

### Shelf life

12 months

**For Cell Culture, further Manufacturing and Research use only. Not approved for human or veterinary use, for application to humans or animals, or for use in therapeutic, clinical or *in vitro* procedures.**

### Support

Life Science Production is a division of Life Science Group Ltd.  
Life Science Production is [ISIA Traceability Certified](#)  
Life Science Group Ltd is an ISO 9001:2015 Certified company

To learn more, contact us:  
Telephone: +44 (0) 1234 889180  
Email: [sales@lifesciencegroup.co.uk](mailto:sales@lifesciencegroup.co.uk)

## QUALITY - TRACEABILITY - INTEGRITY

### Life Science Production

Unit 14, Crowhill Farm, Ravensden Road, Wilden, Bedford, Bedfordshire MK44 2QS  
Tel: 01234 889180 Email: [sales@lifesciencegroup.co.uk](mailto:sales@lifesciencegroup.co.uk) Web: [www.lifesciencegroup.co.uk](http://www.lifesciencegroup.co.uk)

Life Science Production, Antibody Production Services and Diagnostic Kit Services are divisions of Life Science Group Ltd.

Company Registered in England Company Number 06606016

